


Quick Guide EMMA™

World's first self contained Emergency Capnometer

1. EMMA works with all patient connections





2. EMMA ensures proper ventilation & intubation in any environment




3. EMMA is easy to read


Correct Ventilation/Intubation



Warning Hypoventilation





Warning Hyperventilation





Warning incorrect intubation

4. EMMA for easy CPR follow-up



Correct output and RR during CPR



ROSC Return of Spontaneous Circulation

5. EMMA – plug-in and measure



"An EMMA in every emergency bag!"



PHASEIN AB, Svärdvägen 15, SE-182 33 Danderyd, Sweden
 phone: +46 8 544 98 150, fax: +46 8 544 98 169
 info@phasein.com - www.phasein.com